



Date: 28<sup>th</sup> January, 2022

**TO WHOMSOEVER IT MAY CONCERN**

We, **Meril Diagnostics Pvt. Ltd.**, located at Second Floor, D1-D3, Meril Park, Survey No. 135/2/B, 174/2, Muktanand Marg, Chala, Vapi - 396191, Gujarat, India, manufacturer of MERISCREEN COVID-19 Antigen Detection Test declares that our product MERISCREEN COVID-19 Antigen Detection Test detects below mentioned variants of concern of SARS-CoV-2.

- **B.1.1.7(Alpha)**
- **B.1.351(Beta)**
- **P.1(Gamma)**
- **B.1.617.2(Delta)**
- **C.37(Lambda)**
- **B.1.1.529(Omicron)**

Sincerely,



Ram Kanoje  
HEAD - Quality Assurance  
Meril Diagnostics Pvt. Ltd.